








## DIAGNOSTIC



**CS Type**

**SS Type**

## Mouth Mirrors & Handles

Mouth Mirror Back Surface Plane (CS Type)		#3 Ø 20 mm <b>301/3 cs</b> #4 Ø 22 mm <b>301/4 cs</b> #5 Ø 24 mm <b>301/5 cs</b>
Mouth Mirror Back Surface Plane (SS Type)		#3 Ø 20 mm <b>301/3 ss</b> #4 Ø 22 mm <b>301/4 ss</b> #5 Ø 24 mm <b>301/5 ss</b>
Mouth Mirror Front Surface Plane (CS Type)		#3 Ø 20 mm <b>302/3 cs</b> #4 Ø 22 mm <b>302/4 cs</b> #5 Ø 24 mm <b>302/5 cs</b>
Mouth Mirror Front Surface Plane (SS Type)		#3 Ø 20 mm <b>302/3 ss</b> #4 Ø 22 mm <b>302/4 ss</b> #5 Ø 24 mm <b>302/5 ss</b>
Mouth Mirror Back Surface Magnifying (CS Type)		#3 Ø 20 mm <b>303/3 cs</b> #4 Ø 22 mm <b>303/4 cs</b> #5 Ø 24 mm <b>303/5 cs</b>
Mouth Mirror Back Surface Magnifying (SS Type)		#3 Ø 20 mm <b>303/3 ss</b> #4 Ø 22 mm <b>303/4 ss</b> #5 Ø 24 mm <b>303/5 ss</b>
Mouth Mirror Front Surface Double Sided		#3 Ø 20 mm <b>304/3</b> #4 Ø 22 mm <b>304/4</b> #5 Ø 24 mm <b>304/5</b>



**300/1 cs**  
**300/1 ss**  
(mm 120)  
Octagonal handle



**300/2 cs**  
**300/2 ss**  
(mm 120)  
ø mm 6



**300/3 cs**  
**300/3 ss**  
(mm 120)



**300/4 cs**  
**300/4 ss**  
(mm 125)  
Hollow handle  
ø mm 8



**300/5 cs**  
**300/5 ss**  
(mm 125)  
Hollow handle



**300/6 cs**  
**300/6 ss**  
(mm 130)  
Hollow handle

Coton & Dressing Tweezers

DIAGNOSTIC



**01-5001**(mm 130)  
**01-5002**(mm 150)  
Perry



**01-5003** (mm 150)  
Flagg



**01-5004** (mm 160)  
Meriam



**01-5005** (mm 150)  
College



**01-5006** (mm 160)  
Meriam  
Diamond coated



**01-5007**  
(mm 130)  
Perry



**01-5008**  
(mm 150)  
Flagg



**01-5009**  
(mm 160)  
Meriam



**01-5010**  
(mm 150)  
College



**01-5011**  
(mm 150)  
College  
Smooth Jaw





**01-5012**  
(mm 150)  
College  
Light pattern



**01-5013**  
(mm 150)  
College  
With lock

# DIAGNOSTIC

## Explorers Single End




**355/2 S2** Handle Type 2  
**355/2 S1** Handle Type 1

## Explorers Double End



**360/1R** Handle Type 1  
**360/1** Handle Type 2




3-6-9-12

12-9-6-3

## Periodontal Probes

## DIAGNOSTIC



**348/1 S1**  
**348/1 S2**  
**348/1 H2**  
CP12S-O.M.S.

**348/2 S1**  
**348/2 S2**  
**348/2 H2**  
CP11

**348/3 S1**  
**348/3 S2**  
**348/3 H2**  
CP12

**348/4 S1**  
**348/4 S2**  
**348/4 H2**  
CP15

**348/5 S1**  
**348/5 S2**  
**348/5 H2**  
WHO

**348/6 S1**  
**348/6 S2**  
**348/6 H2**  
CP10




**348/1 S2**  
6 mm  
Handle  
Type 2

**348/1 S1**  
6 mm  
Handle  
Type 1



**348/1 H2**  
Hollow Handle 8 mm



**346/1 S1**  
**346/1 S2**  
**346/1 H2**  
Williams

**346/2 S1**  
**346/2 S2**  
**346/2 H2**  
Merritt A

**346/3 S1**  
**346/3 S2**  
**346/3 H2**  
Merritt E

**346/4 S1**  
**346/4 S2**  
**346/4 H2**  
Goldman Fox


**346/5 S1**  
**346/5 S2**  
**346/5 H2**  
University of Michigan

**346/6 S1**  
**346/6 S2**  
**346/6 H2**  
26 G

**346/7 S1**  
**346/7 S2**  
**346/7 H2**  
Mercury

# DIAGNOSTIC

## Probe Double End




**367/3 S1**  
**367/3 S2**  
Williams



**367/16 S1**  
**367/16 S2**  
Williams



**367/23 S1**  
**367/23 S2**  
Williams



**367/17 S1**  
**367/17 S2**  
Williams










**365/3 S1**  
**365/3 S2**  
University  
of Michigan



**365/16 S1**  
**365/16 S2**  
University  
of Michigan




**365/23 S1**  
**365/23 S2**  
University  
of Michigan



**365/23 S2**  
6 mm  
Handle  
Type 2




**365/23 S1**  
6 mm  
Handle  
Type 1




**366/16 S1**  
**366/16 S2**  
Goldman-Fox



**366/23 S1**  
**366/23 S2**  
Goldman-Fox



**366/17 S1**  
**366/17 S2**  
Goldman-Fox



**368/1 S1**  
**368/1 S2**  
Goldman-Fox  
Williams



**370/1 S1**  
**370/1 S2**  
CP11/23



**370/2 S1**  
**370/2 S2**  
CP12/23



**370/3 S1**  
**370/3 S2**  
CP15/23

